


Don't delay.

Use of this injection will not



Paediatric Emergency
Hydrocortisone
Intramuscular Injection Kit

Contents:

Plastic bag to use as surface;

2 x Hydrocortisone 100mg/mL ampoule;

2 x Syringes; 2 x Blue needles;





Scan QR code to access audio and video recording of leaflet

Tap top of ampoule so all fluid is at the bottom.



Hold ampoule in one hand, with thumb of other hand on spot on the neck.

Snap neck with thumb.

Attach needle to syringe and remove cap.





Place syringe in ampoule with needle tip touching the bottom.





Pull syringe plunger towards you to collect the required volume of fluid:




Child less than 1 year = 0.25mL (25mg)

1-5 years = 0.5mL (50mg) Over 5 years = 1mL (100mg)









Hold syringe with needle up and tap to let air bubbles rise.

Push plunger to remove air.




Locate injection site in middle outer thigh.







Call 999 and tell them emergency adrenal crisis.



Dispose of kit container in sharps bin in ambulance, hospital, or pharmacy.

Contact details: